

## Question-Q13- RespiratoryRateAndAudioBasedTracking

**Question: Can Respiratory Rate + Movement accurately track sleep stages, and will a pillow-located sensor increase the accuracy of this type of monitoring?**

---

### Answers

#### TLDR

### Respiratory Rate and Pillow-Based Sleep Monitoring: Executive Summary

#### Research Question

Can respiratory rate + movement accurately track sleep stages, and will a pillow-located sensor increase the accuracy of this type of monitoring?

#### Key Findings

##### Respiratory Rate as Sleep Stage Indicator

- **Scientific Validity:** Research confirms respiratory patterns reliably change across sleep stages
  - NREM sleep: Regular, decreased respiratory rate
  - REM sleep: Faster breathing with higher variability
  - Respiratory rate variability (RRV) correlates with sleep depth (lowest during deep sleep)
- **Accuracy:** When combined with movement data, respiratory monitoring shows up to 87% agreement with polysomnography (the clinical gold standard)
- **Respiratory Patterns:** 97.5% of breathing sound events occur during inspiration, making inspiratory sounds particularly valuable for analysis

# Pillow-Based Sleep Tracking Technologies

Several companies and research teams are developing promising solutions:

1. **The Sleep Company's SensAI Sleep Analyzer Pillow**
    - Non-intrusive sensor strip placed under pillow and mattress
    - Tracks sleep duration, quality, body movements, heart rate, and breathing rate
    - Uses low-energy radar technology to avoid wearable discomfort
  2. **Smart Pillow with FB-TENG Technology**
    - Incorporates flexible and breathable triboelectric nanogenerator arrays
    - Generates electrical signals based on head movement and pressure
    - Self-powering system for real-time tracking
  3. **Temperature and Humidity Sensing Pillow**
    - Strategic sensor placement within the pillow
    - Detects body temperature and night sweats
    - Bluetooth transmission to cloud-based health sensing system
  4. **Pressure and Sound Sensor Integration**
    - Systems using Raspberry Pi with both pressure and sound sensors
    - Pressure sensors track sleep duration and REM cycles
    - Sound sensors detect sleep-disordered breathing and snoring
  5. **AJProTech's Smart Pillow**
    - Multi-sensor approach including microphones and gyroscopes
    - LED light technology for sleep induction and waking
    - White noise features and Bluetooth connectivity
- 

## Advantages of Pillow-Based Monitoring

Pillow-based technologies offer significant benefits over smartphone/wearable approaches:

1. **Consistent Sensor Placement:** Fixed position eliminates placement variability
2. **Environmental Noise Reduction:** Physical filtering of ambient sounds
3. **User Isolation:** Better separation of individual data in shared beds
4. **Multi-Modal Sensing:** Combines acoustic, pressure, movement, and temperature data
5. **User Comfort:** Unobtrusive, non-wearable solution enhances compliance

## Conclusion

Respiratory rate monitoring combined with movement detection shows strong potential for accurate, non-invasive sleep tracking. Pillow-based systems represent a promising approach by addressing key limitations of current technologies while maintaining user comfort. While several research prototypes and commercial products are exploring this space, significant innovation opportunities remain for developing comprehensive sleep monitoring solutions.

## Market Status

As of March 2025, while various pillow-based prototypes exist, fully integrated commercial systems with microphones and comprehensive sensor arrays are still emerging, suggesting this remains an area with substantial growth potential.

Respiratory rate is indeed a powerful indicator of sleep stages, and several technologies are leveraging this insight through pillows with integrated sensors. Here's what the research shows:

## Respiratory Rate as a Sleep Stage Indicator

The relationship between respiratory patterns and sleep stages is well-established in research:

- Breathing rate becomes regular and decreases upon entering NREM sleep, while REM sleep is associated with faster breathing rate and higher variability<sup>2</sup>
- Respiratory rate variability (RRV) is greatest during wakefulness and correlates with sleep depth, being lowest during N3 (deep sleep)<sup>6</sup>
- The CardioRespiratory Sleep Staging (CReSS) algorithm demonstrates that combining heart rate variability with respiration data can effectively identify sleep stages (wake, light sleep, deep sleep, and REM)<sup>1</sup>
- A sound-based sleep staging system achieved 87% agreement with polysomnography (the gold standard) by analyzing breathing and movement sounds<sup>2</sup>

Most breathing sound events (97.5%) occur during inspiration, making inspiratory sounds particularly valuable for analysis<sup>7</sup>. Snoring patterns also vary with sleep stages - snoring index is higher in slow wave sleep compared to S2 and REM sleep<sup>7</sup>.

## Pillow-Based Sleep Tracking Technologies

Several companies are developing pillow-based sleep tracking systems that address the limitations you mentioned:

### The Sleep Company's SensAI Sleep Analyzer Pillow

- Uses a non-intrusive sensor strip placed under the pillow and mattress
- Tracks sleep duration, quality (including deep and REM stages), body movements, heart rate, and breathing rate

- Employs low-energy radar technology to avoid wearable discomfort
- Data is transmitted to an interactive app for detailed insights<sup>5</sup>

### **Smart Pillow with FB-TENG Technology**

- Developed by researchers Ding Li and Zhong Lin Wang
- Incorporates flexible and breathable triboelectric nanogenerator arrays on the pillow surface
- Generates electrical signals based on head movement and pressure distribution
- Self-powering system that provides accurate, real-time tracking of head position<sup>810</sup>

### **Temperature and Humidity Sensing Pillow**

- Implants temperature and humidity sensors at strategic positions within the pillow
- Infers body temperature and detects night sweats
- Data is transmitted via Bluetooth to a cloud-based health sensing system<sup>4</sup>

### **AJProTech's Smart Pillow**

- Uses a combination of sensors to track sleep
- Incorporates LED light technology (red light for sleep induction, gradually brightening light for waking)
- Includes white noise features and Bluetooth connectivity<sup>9</sup>

These pillow-based technologies offer significant advantages over smartphone microphone-based systems by:

1. Providing consistent sensor placement
2. Reducing interference from environmental sounds
3. Better isolating individual user data in shared beds
4. Combining multiple sensing modalities (sound, pressure, movement, temperature)

The evolution of these technologies suggests that your insight about using the pillow as a hub for sleep sensors is already being pursued by researchers and companies with promising results.

### **Citations:**

1. <https://pmc.ncbi.nlm.nih.gov/articles/PMC8314617/>
2. <https://pmc.ncbi.nlm.nih.gov/articles/PMC6128888/>
3. <https://pillow.app/article/snore-detection-and-tracking>
4. <https://pmc.ncbi.nlm.nih.gov/articles/PMC6263409/>

5. <https://www.expresscomputer.in/news/the-sleep-company-bets-big-on-ai-powered-comfort-launches-worlds-first-smartgrid-sleep-analyser-pillow-sensai/119642/>
6. [https://hsrc.himmelfarb.gwu.edu/context/smhs\\_surgery\\_facpubs/article/1418/viewcontent/Respiratory\\_rate\\_variability\\_in\\_sleeping\\_adults\\_without\\_obstructive\\_sleep\\_apnea.pdf](https://hsrc.himmelfarb.gwu.edu/context/smhs_surgery_facpubs/article/1418/viewcontent/Respiratory_rate_variability_in_sleeping_adults_without_obstructive_sleep_apnea.pdf)
7. <https://jcsn.aasm.org/doi/10.5664/jcsn.5588>
8. <https://www.biotechniques.com/neuroscience/comfortable-and-clever-the-smart-pillow-that-tracks-your-sleep-behavior/>
9. <https://ajprotech.com/selected-work/health-wellness/sleep-tracking-pillow.html>
10. <https://www.acs.org/pressroom/newsreleases/2022/june/tracking-sleep-with-a-self-powering-smart-pillow.html>
11. <https://www.ncbi.nlm.nih.gov/books/NBK526132/>
12. <https://ouraring.com/blog/breathing-rate-while-sleeping/>
13. <https://www.sleepfoundation.org/sleep-apnea/sleep-respiratory-rate>
14. <https://ieeexplore.ieee.org/document/8512446/>
15. <https://pmc.ncbi.nlm.nih.gov/articles/PMC11041471/>
16. <https://jamanetwork.com/journals/jamaotolaryngology/fullarticle/2791067>
17. <https://www.whoop.com/us/en/thelocker/understanding-respiratory-rate-what-it-is-whats-normal-and-why-you-should>
18. <http://ieeexplore.ieee.org/document/8688505/>
19. <https://www.mayoclinic.org/tests-procedures/polysomnography/about/pac-20394877>
20. <https://www.sleepfoundation.org/stages-of-sleep>
21. [https://www.atsjournals.org/doi/pdf/10.1164/ajrccm-conference.2016.193.1\\_MeetingAbstracts.A5983?download=true](https://www.atsjournals.org/doi/pdf/10.1164/ajrccm-conference.2016.193.1_MeetingAbstracts.A5983?download=true)
22. <https://www.sleepapnea.org/sleep-health/sleep-respiratory-rate/>
23. <https://pmc.ncbi.nlm.nih.gov/articles/PMC10654909/>
24. [https://intercom.help/pillow\\_app/en/articles/3733515-pillow-didn-t-record-any-sounds-during-sleep](https://intercom.help/pillow_app/en/articles/3733515-pillow-didn-t-record-any-sounds-during-sleep)
25. <https://www.teslasmart.com/tesla-smart-pillow>
26. <https://newsroom.mattressfirm.com/press/mattress-firm-empowers-consumers-to-improve-well-being-with-new-sleep-tracking-app/>
27. <https://pillow.app/pillow-v4-5-audio-recordings-and-machine-learning>
28. <https://pillow.app>
29. <https://vertu.com/lifestyle/top-sleep-tracking-innovations-2025/>
30. [https://www.reddit.com/r/AppleWatch/comments/a6zx8n/sleep\\_apps\\_microphone\\_access\\_paranoid/](https://www.reddit.com/r/AppleWatch/comments/a6zx8n/sleep_apps_microphone_access_paranoid/)
31. <https://us.sleepace.com/pages/smartpillow?id=938004>
32. <https://www.nytimes.com/wirecutter/reviews/best-sleep-trackers/>

33. [https://intercom.help/pillow\\_app/en/articles/3733541-why-does-pillow-need-access-to-the-microphone](https://intercom.help/pillow_app/en/articles/3733541-why-does-pillow-need-access-to-the-microphone)
  34. <http://motionpillow.com>
  35. <https://www.nature.com/articles/s41746-021-00493-6>
  36. <https://pmc.ncbi.nlm.nih.gov/articles/PMC9996414/>
  37. <https://sleepdoctor.com/sleep-apnea/respiratory-rate-while-sleeping>
  38. <https://aasm.org/simple-analysis-of-breathing-sounds-while-awake-can-detect-obstructive-sleep-apnea/>
  39. <https://apps.apple.com/gb/app/pillow-sleep-tracker/id878691772>
  40. <https://apps.apple.com/by/app/pillow-sleep-tracker/id878691772>
- 

Answer from Perplexity: [pplx.ai/share](https://pplx.ai/share)

---

## Key Points

- Research suggests respiratory rate is a key indicator for sleep stages, especially when combined with movement data.
  - It seems likely that accuracy issues from environmental noise can be mitigated with a microphone under the pillow, potentially including accelerometers or pressure sensors.
  - The evidence leans toward some companies exploring pillow-based sleep trackers, but dedicated devices with microphones are not widely available yet.
- 

## Respiratory Rate and Sleep Stages

Respiratory rate changes during different sleep stages, with slower, more regular breathing in deep sleep and irregular, faster breathing in REM sleep. While respiratory rate alone can provide insights, combining it with movement detection enhances accuracy for measuring sleep stages. This approach is used in many consumer sleep trackers, like smartphone apps, to estimate sleep quality and stages.

## Microphone Under the Pillow: Potential and Current Developments

Placing a microphone under the pillow could reduce environmental noise interference, improving accuracy for sleep tracking. Adding accelerometers or pressure sensors could further detect movement, offering a comprehensive sleep analysis. While no widely available dedicated devices with microphones in pillows exist, some companies, like TESLA with their Smart Pillow, are developing sleep monitoring features, though not specifically mentioning microphones. Research also shows ongoing innovation, such as self-powering smart pillows tracking head movements, but these focus more on motion than sound.

## Conclusion

Your observations highlight promising avenues for sleep tracking technology, particularly the integration of respiratory rate and movement data with innovative placement like under the pillow. While current developments are promising, more dedicated devices may emerge in the future.

---

## Survey Note: Detailed Analysis of Sleep Tracking Technologies

This section provides a comprehensive exploration of sleep tracking technologies, focusing on the role of respiratory rate in determining sleep stages and the potential for pillow-based systems with microphones, accelerometers, and pressure sensors. The analysis is grounded in recent research and product developments, aiming to offer a thorough understanding for researchers and enthusiasts alike.

### Respiratory Rate as a Sleep Stage Indicator

Sleep is characterized by distinct stages—awake, light sleep, deep sleep, and REM sleep—each with unique physiological signatures, including breathing patterns. Research suggests that respiratory rate is a powerful indicator for distinguishing these stages. For instance, during deep sleep, breathing tends to be slower and more regular, as noted in a 2023 article from the Sleep Foundation ([Understanding Your Respiratory Rate While Sleeping](#)). In contrast, REM sleep often features irregular and faster breathing, as highlighted in a 2024 SleepApnea.org article ([Sleep Respiratory Rate](#)).

A 2016 study published in the National Center for Biotechnology Information (NCBI) examined respiratory rate variability (RRV) in sleeping adults without obstructive sleep apnea, finding greater variability during wakefulness compared to sleep stages, including REM ([Respiratory rate variability in sleeping adults without obstructive sleep apnea](#)). This suggests that respiratory rate can help differentiate sleep from wakefulness, but accuracy for specific stages may require additional data.

The combination of respiratory rate with movement detection appears to enhance accuracy. For example, WHOOP, a fitness tracker company, uses minute-by-minute alterations in respiratory rate within its sleep staging algorithm, as mentioned in a 2021 article ([Respiratory Rate: What's Normal & Why You Should Track It](#)). This approach aligns with the user's observation that the app likely relies on breathing sounds and movement detection, indicating a synergistic effect. A 1982 study in Thorax further supports this, showing that ventilation falls during sleep, with the greatest reduction in REM, and breathing patterns become more rapid and shallow ([Respiration during sleep in normal man](#)).

However, respiratory rate alone may not be sufficient for precise sleep stage classification, as traditional polysomnography (PSG) uses EEG, EMG, and EOG alongside respiratory and cardiac measures. The evidence leans toward combining respiratory rate with movement for better results, especially in consumer devices where simplicity is key.

## **Mitigating Accuracy Issues with Pillow-Based Technology**

The user's insight into placing a microphone under the pillow to mitigate accuracy issues from environmental noise and multiple sleepers is particularly noteworthy. This approach could isolate breathing sounds more effectively, reducing interference from external sources. Adding an accelerometer or pressure sensor under the pillow could further detect movement, providing a dual-modality system for comprehensive sleep analysis.

Current sleep tracking apps, such as Pillow and Sleep Cycle, often use the smartphone's microphone when placed near the bed or under the pillow. For instance, Sleep Cycle's patented sound analysis technology uses the microphone to detect bed sheet movements, as detailed in a 2018 article ([Sound analysis sleep tracking - Patented - Sleep Cycle alarm clock](#)). However, these are not dedicated devices, and placing a smartphone under a pillow can lead to overheating, as noted in Sleep Cycle's support documentation ([Phone Placement with the Sleep Cycle App](#)).

Dedicated devices are less common, but some innovations are emerging. The TESLA Smart Pillow, as described on their website, includes sleep monitoring among its features, with data viewable via a mobile app ([TESLA Smart Pillow](#)). However, it does not explicitly mention a microphone, suggesting it may rely on other sensors like temperature or movement. A 2022 study in the American Chemical Society discussed a self-powering smart pillow using triboelectric nanogenerators to track head movements, focusing on motion rather than sound ([Tracking sleep with a self-powering smart pillow](#)).

An interesting development is mentioned in a case study by AJProTech, which describes a "Smart Pillow" with a microphone to analyze breathing patterns and detect snoring, as well as gyroscopes and accelerometers for movement detection ([Case study | Sleep Tracking Pillow |](#)



[IoT Development](#)). However, the specific page was inaccessible, indicating potential limitations in publicly available information.

Research also points to ongoing efforts in this area. A 2019 PMC article discussed consumer sleep apps, noting that phone-near-pillow apps struggle to distinguish sleep stages but can reasonably detect sleep versus wake, similar to wrist actigraphy ([Consumer Sleep Apps: When it Comes to the Big Picture, it's All About the Frame](#)). This suggests that a dedicated pillow-based system with a microphone could improve stage differentiation, especially if combined with movement sensors.

## Companies and Future Prospects

While no widely available commercial products feature a microphone embedded in a pillow for sleep tracking, several companies and research groups are exploring related technologies. The TESLA Smart Pillow is one example, though its sleep monitoring method is unclear. Another is the Sleepace Sleep Dot, mentioned in a 2017 Mashable article, which sticks to the pillowcase to track movement and sleep cycles, but does not specify a microphone ([From vibrating pillowcases to smart pajama belts, sleep tech is flooding the market](#)).

Chinese researchers, as reported in I'MNOVATION, have developed a smart pillow using triboelectric nanogenerators to track head movements, offering high precision but focusing on motion rather than sound ([Smart Pillows: Sweet Dreams Are Made of This](#)). This indicates a trend toward integrating advanced sensors into pillows, though microphone-based systems are less documented.

The concept of a pillow as a hub for sleep sensors is gaining traction, with potential for future devices to include microphones, accelerometers, and pressure sensors. For instance, Amazon filed a patent in 2024 for a smart pillow to combat sleep apnea, suggesting interest in pillow-based solutions, though details on sensors were not specified ([Amazon could be making its own smart pillow to combat sleep apnea – but I'm not convinced](#)).

Given these developments, it seems likely that more companies will explore pillow-based sleep trackers with microphones, especially as the demand for non-invasive, comfortable sleep monitoring grows. However, as of March 13, 2025, such devices are not yet mainstream, and researchers may need to monitor CES and other tech expos for updates, as seen in recent coverage ([The best sleep tech from CES 2025 so far](#)).

## Conclusion and Implications

The user's observations align with current research, highlighting respiratory rate as a key metric for sleep stage analysis, particularly when combined with movement data. The idea of a microphone under the pillow to reduce noise interference is promising and could lead to more accurate sleep tracking, especially with additional sensors like accelerometers and pressure

detectors. While dedicated devices are not yet widely available, ongoing innovations suggest a future where pillows serve as hubs for comprehensive sleep monitoring, potentially revolutionizing consumer sleep technology.

This analysis underscores the need for further research into integrating multiple sensor modalities in pillow-based systems, ensuring both accuracy and comfort for users. Researchers interested in this field should explore patent databases and scientific literature for emerging technologies, as well as engage with companies like TESLA and Sleepace for potential collaborations.

---

## Key Citations

- [Understanding Your Respiratory Rate While Sleeping](#)
  - [Sleep Respiratory Rate](#)
  - [Respiratory rate variability in sleeping adults without obstructive sleep apnea](#)
  - [Respiration during sleep in normal man](#)
  - [Respiratory Rate: What's Normal & Why You Should Track It](#)
  - [Pillow: Sleep Tracker, Monitor & Alarm Clock](#)
  - [Phone Placement with the Sleep Cycle App](#)
  - [TESLA Smart Pillow](#)
  - [Tracking sleep with a self-powering smart pillow](#)
  - [Sound analysis sleep tracking - Patented - Sleep Cycle alarm clock](#)
  - [Case study | Sleep Tracking Pillow | IoT Development](#)
  - [Consumer Sleep Apps: When it Comes to the Big Picture, it's All About the Frame](#)
  - [From vibrating pillowcases to smart pajama belts, sleep tech is flooding the market](#)
  - [Amazon could be making its own smart pillow to combat sleep apnea – but I'm not convinced](#)
  - [Smart Pillows: Sweet Dreams Are Made of This](#)
  - [The best sleep tech from CES 2025 so far](#)
- 

## Respiratory Rate and Pillow-Based Technologies: Advanced Approaches to Sleep Stage Monitoring

The exploration of respiratory patterns during sleep reveals significant potential for non-invasive sleep stage detection, while emerging pillow-based technologies offer promising solutions for

convenient sleep monitoring. These advancements could revolutionize how we track and understand sleep quality, particularly for individuals seeking alternatives to wearable devices or clinical sleep studies.

## The Relationship Between Respiratory Rate and Sleep Stages

Respiratory patterns undergo distinct changes across different sleep stages, making breathing a valuable biomarker for sleep monitoring. During the transition from wakefulness to non-rapid eye movement (NREM) sleep, there is a progressive reduction in central respiratory drive as sleep deepens from stage N1 through N3<sup>5</sup>. This physiological shift results in measurable changes in breathing patterns that correlate strongly with sleep stage progression.

Research has demonstrated that respiratory rate variability (RRV) – the variation in breathing patterns – decreases consistently with sleep depth. RRV is greatest during wakefulness, decreases during REM sleep, and progressively diminishes through the NREM stages, reaching its lowest point during N3 (deep sleep)<sup>1</sup>. This systematic change provides a reliable indicator for differentiating between sleep stages. Additionally, a person's baseline respiratory rate naturally decreases during sleep, with the slowest rates occurring during deep sleep, while REM sleep is characterized by increased respiratory rate and irregularity<sup>6</sup>.

These respiratory changes during sleep have physiological underpinnings. During NREM sleep, respiration becomes more regular and predominantly falls under metabolic control, whereas during REM sleep, respiration becomes irregular and is influenced more by behavioral factors<sup>5</sup>. The PaCO<sub>2</sub> levels typically increase by 3-7 mmHg during sleep, further reflecting these respiratory adjustments<sup>5</sup>. For healthy adults between 18-65 years, normal respiratory rates range from 12-20 breaths per minute while awake, with rates declining during sleep<sup>6</sup>.

Studies investigating breathing irregularities during REM sleep have found fascinating connections between respiratory patterns and dream content. Research with narcolepsy patients capable of lucid dreaming showed congruence between dream reports and observed respiratory behaviors, suggesting that breathing irregularities during REM sleep may have a cortical or subcortical origin reflecting dream content<sup>10</sup>. This adds another dimension to understanding the relationship between breathing patterns and sleep states.

## Pillow-Based Sleep Tracking Technologies

The limitations of smartphone microphone-based sleep trackers – including variability in placement, environmental sounds, and multiple sleepers – have spurred innovation in pillow-based sleep monitoring. Several promising technologies are emerging in this space, offering potential solutions that address these challenges while maintaining user comfort.

# Pressure and Sound Sensor Integration

Researchers have developed pillow-based sleep tracking devices using Raspberry Pi that incorporate both pressure sensors and sound sensors<sup>3</sup>. This integration allows for monitoring of sleep duration, sleep stages, and detection of snoring. The pressure sensor, connected to a pressure amplifier bag placed under the pillow, targets the detection of sleep duration and REM cycles, while the sound sensor captures sleep-disordered respiration. This approach offers a compelling alternative to wearable devices by utilizing an item that is naturally and involuntarily used during sleep<sup>3</sup>.

# Temperature and Humidity Sensing Systems

Another innovative approach involves implanting temperature and humidity sensors inside the pillow in strategic positions<sup>4</sup>. This smart pillow system identifies the roles of different sensors (main, auxiliary, or environmental temperature) based on value differences, allowing for pattern extraction through statistical analysis. Body temperature can be inferred through a specially designed Fuzzy Logic System when the head position remains stable. The system can also detect night sweats using humidity sensors. Data is transmitted via Bluetooth Low Energy and uploaded to a cloud server for analysis, creating a comprehensive health sensing system<sup>4</sup>.

# Triboelectric Nanogenerator Technology

Advanced smart pillows using flexible and breathable triboelectric nanogenerator (FB-TENG) technology represent another promising development<sup>8</sup>. This technology converts mechanical energy into electrical energy, serving as the basis for pressure sensors. The smart pillow incorporates a layer of FB-TENGs on its surface, forming a porous array of self-powering sensors. This design allows for comfortable recording of head movements during sleep, offering an attractive alternative to medical sleep tests and smartphone trackers<sup>8</sup>.

# Combined Acoustic and Motion Detection Systems

While not specifically pillow-based, some integrated systems demonstrate the effectiveness of combining acoustic and motion detection for sleep monitoring. One such system uses a special electret microphone and an inertial measurement unit (IMU) to detect heartbeats, breathing, snoring, sleeping positions, and movements<sup>2</sup>. The system analyzes sleep using acoustic information from the microphone to identify breathing events and heartbeats, while the IMU provides data on patient activity and positions. The information from both sensors is then fused to detect sleep events, showing promise for improved sleep monitoring<sup>2</sup>.

# Commercial Applications and Consumer Technologies

In the commercial space, several technologies have emerged that apply these principles. The Nanit Breathing Band, although designed for infants, demonstrates how pattern recognition can be used for respiratory monitoring without direct sensors<sup>7</sup>. The band works with a baby monitor to detect breathing motion using AI-powered technology, tracking breaths per minute without sensors, electronics, or radars directly on the child<sup>7</sup>.

The "Pillow" app for Apple devices represents a software-based approach that analyzes sleep cycles using data from Apple Watch, iPhone, or iPad<sup>9</sup>. While not a physical pillow solution, it demonstrates the commercial viability of sleep tracking technology that incorporates respiratory and movement data.

## Implications for Future Technology Development

The convergence of respiratory monitoring and pillow-based technologies offers significant potential for advancing sleep tracking. By placing microphones under pillows, many of the limitations of smartphone-based sleep trackers could be mitigated. The pillow provides a natural fixed position for the microphone, reducing placement variability, while also offering some filtering of environmental sounds. Furthermore, the proximity to the sleeper's head creates optimal conditions for detecting breathing sounds.

Combining this acoustic approach with pressure sensors or accelerometers in the pillow could yield a comprehensive sleep monitoring solution. The pressure sensor would detect movement and potentially heart rate, while the microphone would capture respiratory patterns. Given the research showing how respiratory rate and its variability change across sleep stages, this combined approach could provide reliable sleep stage classification without the need for electroencephalography (EEG)<sup>156</sup>.

For individuals sharing a bed, the localized nature of pillow-based monitoring would help isolate the target sleeper's signals from those of bed partners, addressing another major limitation of smartphone-based solutions. Additionally, the unobtrusive nature of pillow-based technology eliminates the discomfort associated with wearable devices, potentially improving user adherence.

## Conclusion

Respiratory rate and its variability represent powerful indicators of sleep stage, with distinct patterns emerging across the wake-sleep cycle and different sleep stages. This physiological relationship forms a solid foundation for non-invasive sleep tracking technologies. Pillow-based monitoring systems offer a particularly promising approach, combining user comfort with optimal positioning for detecting respiratory and movement signals.

While several research prototypes and commercial products are already exploring this space, there remains significant room for innovation. Future developments might include machine

learning algorithms that become increasingly accurate at distinguishing sleep stages based on respiratory patterns, or smart pillows that adjust their properties based on detected sleep states to optimize sleep quality. As these technologies continue to evolve, they promise to make sophisticated sleep tracking more accessible, comfortable, and accurate for everyday users.

## Citations:

1. <https://pmc.ncbi.nlm.nih.gov/articles/PMC5027356/>
2. <https://pubmed.ncbi.nlm.nih.gov/25571592/>
3. <https://scholarworks.lib.csusb.edu/cgi/viewcontent.cgi?article=2790&context=etd>
4. <https://pmc.ncbi.nlm.nih.gov/articles/PMC6263409/>
5. <https://pmc.ncbi.nlm.nih.gov/articles/PMC2876696/>
6. <https://www.sleepapnea.org/sleep-health/sleep-respiratory-rate/>
7. <https://www.nanit.com/products/nanit-breathing-band>
8. <https://www.biotechniques.com/neuroscience/comfortable-and-clever-the-smart-pillow-that-tracks-your-sleep-behavior/>
9. <https://apps.apple.com/us/app/pillow-sleep-tracker/id878691772>
10. <https://www.nature.com/articles/s41598-018-21067-9>
11. <https://www.nature.com/articles/1300146>
12. <https://jcsn.aasm.org/doi/10.5664/jcsn.6778>
13. <https://www.frontiersin.org/journals/sleep/articles/10.3389/frsle.2023.1145620/full>
14. <https://www.withings.com/eu/en/sleep-analyzer>
15. <https://www.babysensemonitors.com/products/babysense-7-breathing-monitor>
16. <https://www.nytimes.com/wirecutter/reviews/best-sleep-trackers/>
17. <https://pmc.ncbi.nlm.nih.gov/articles/PMC459437/>
18. <https://fetalplus.com/products/fetalplus-baby-sleep-breathing-monitor-wireless-smart-gateway>
19. <https://www.advancedsciencenews.com/a-smart-pillowcase-that-tracks-sleep/>
20. [https://www.reddit.com/r/AppleWatch/comments/l1idn1/how\\_to\\_use\\_pillow\\_for\\_sleep\\_tracking/](https://www.reddit.com/r/AppleWatch/comments/l1idn1/how_to_use_pillow_for_sleep_tracking/)
21. <https://www.nationaljewish.org/njh/media/pdf/pdf-vipin-malik-publication-respiratory.pdf>
22. <https://shop.sense-u.com/en-na/products/sense-u-baby-monitor-pro>
23. <https://us.sleepace.com/pages/smartpillow?id=938004>
24. <https://www.cnet.com/health/sleep/best-sleep-tracker/>
25. <https://www.teslasmart.com/tesla-smart-pillow>
26. <https://www.eightsleep.com>
27. <https://www.nytimes.com/wirecutter/blog/somnox-sleep-robot-review/>
28. [https://www.reddit.com/r/adhdwomen/comments/11vm930/underpillow\\_speakers/](https://www.reddit.com/r/adhdwomen/comments/11vm930/underpillow_speakers/)

29. <https://pubmed.ncbi.nlm.nih.gov/35800029/>
  30. <https://pillow.app>
  31. <https://store.bruker.com/products/respiration-pillow-sensor>
  32. <https://www.mayoclinic.org/tests-procedures/polysomnography/about/pac-20394877>
  33. <https://academic.oup.com/sleep/article/36/7/1069/2453878>
- 

Answer from Perplexity: [pplx.ai/share](https://pplx.ai/share)